

EFFECTIVE JANUARY 1, 2013



Remove the shrink band before filling the bottle

Please fill to "shoulder" of bottle

100 ml line etched on bottle

**UNDERFILLED AND OVERFILLED  
BOTTLES WILL BE REJECTED !!!**

- DON'T TOUCH THE INSIDE OF THE LID or BOTTLE
- DON'T RINSE OUT THE BOTTLE
- KEEP THE SAMPLE COOL WITH ICE OR ICE PACKS AND TRANSPORT TO THE LAB WITHIN 28 HRS.

## CONTACT US

PHONE 832.393.3939

E-MAIL [WaterLab.Info@houstontx.gov](mailto:WaterLab.Info@houstontx.gov)

[WWW.HOUSTONTX.GOV/HEALTH/WATER.HTML](http://WWW.HOUSTONTX.GOV/HEALTH/WATER.HTML)